



# CO 8000 Biowave Personal Cell Density Meter

- Simple
- Reliable
- Portable



## Rapid, low cost cell density measurement, anytime, anywhere

The Biowave CO 8000 from Biochrom is an easy-to-use instrument to measure the density of cells in suspension (e.g. yeast cells and bacterial cultures) by measuring Absorbance (OD) at 600nm

Completely portable, this hand-held instrument has a rechargeable battery, so it can be used in incubation cabinets, under anaerobic conditions, or in fume hoods and safety cabinets.

- LED Light Source for stable results and no maintenance.
- No warm up time.
- Saves last 99 results for download to PC or Printer.
- Accepts 10mm cuvettes, or 16mm tubes (Ehrlenmeyer side-arm flasks).
- Adapters available for 10mm and 12mm tubes.
- Easy to clean.

The CO 8000 gives a lifetime of stable and repeatable results for a very affordable price!



## Technical Information

Light Source	LED
Absorbance Range	-0.3 to 1.99 A
Wavelength Range	600nm
Accuracy	± 0.05 A at 1 A using neutral density filters
Reproducibility	± 0.02 A at 1 A using cuvettes
Bandwidth	40nm
Stray Light	<1%T
Test Tubes	16mm standard round tube. 10mm and 12mm with adapter
Cuvettes	10 x 10mm, semi micro, macro
Display	LCD
Output	RS232 digital
Power	110/240 V via mains transformer/NiMH battery
WxDxH	150mm x 180mm x 60mm, 5.9" x 7.1" x 2.4"
Weight	0.6kg, 1.3lbs
Warranty	1 year

CE marked and complies with relevant legislation, including EMC and low voltage directives.

## Ordering Information

	Order Code
CO 8000 Personal Cell Density Meter mains/rechargeable battery	80-3000-45
Adapter for 10mm and 12mm tubes	80-3000-57
S2000P Serial Printer 40 column incl. cable	80-3000-94
RS232 Serial Cable	80-3001-00

For more information and technical specifications see our websites.  
[www.biochrom.co.uk](http://www.biochrom.co.uk) or [www.biochrom-US.com](http://www.biochrom-US.com)

Distributors worldwide



Biochrom Ltd, 22 Cambridge Science Park, Milton Road, Cambridge, CB4 0FJ England  
Tel: +44 (0)1223 423723 Fax: +44 (0) 1223 420164  
[enquiries@biochrom.co.uk](mailto:enquiries@biochrom.co.uk) [www.biochrom.co.uk](http://www.biochrom.co.uk)

Biochrom-US  
Tel: +1 877-246-2476 Fax: +1 508-429-5732  
[sales@biochrom-us.com](mailto:sales@biochrom-us.com) [www.biochrom-US.com](http://www.biochrom-US.com)